## The KaVo K-POWERgrip: Unbeaten in technology and design.

Ravo EWI

Where absolute top performan-ces are required, there is no alternative to KaVo KPOWERgrip.

The new motor is dynamic and powerful and has exceptional high torque of up to 7 Ncm. This makes the K-POWERgrip extre-mely flexible in use and suitable for all areas of work and materi-

The virtually noiseless and vibra-tion-free induction motor assists concentration. The new, ultra-slim cable and reduced weight gives the user additional free-dom of movement and superb handling.

Particularly important: The new grip surfaces, ensure that the K-POWERgrip sits more comfor-tably and firmly in the hand.

The KaVo K-POWERgrip. The new king among laboratory motors. Its operation is just as good as its looks.



The K-Control, the universal control unit for most KaVo handpieces.





Not only does it look great but also feels great with the superb new grip: the K-POWERgrip in the space-saving handpiece support.



K-Control available with knee, foot and table top control.

Mat. No.: Speed range: Maximum torque: Chuck size Fan: Drive system: Housing shape:

1.002.3589 1.000 - 50.000 min<sup>-1</sup> 7 Ncm 2,35 mm, 3,00 mm, 1/8" ø 0,02 mm without additional cooling brushless DC motor ergonomic shape for the hand with grip inlays





## The KaVo K-POWERgrip. Out in front in every respect.



Always a hundred percent secu-re grip. More freedom of movement More power for every work during work. Situation. re grip.

## 20% cooler handpiece.

The result of improved motor technology is better performance.



Because the K-POWERgrip has a completely new grip surface. The new flexible K-POWERgrip This a new motor technology handpiece cable is ultra-slim, provides more power and torque. ultra-light and flexible.

50



KaVo K-POWERgrip

## Get a grip on real power!







KaVo. Dental Excellence.